The Nelson Mandela AFrican Institution of Science and Technology

https://dspace.mm-aist.ac.tz

Life sciences and Bio-engineering

Research Articles [LISBE]

2021-06-09

TB or not TB? Definitive determination of species within the Mycobacterium tuberculosis complex in unprocessed sputum from adults with presumed multidrug-resistant tuberculosis

Mbelele, Peter

Wiley

https://doi.org/10.1111/tmi.13638 Provided with love from The Nelson Mandela African Institution of Science and Technology VOLUME 26 NO 9 PP 1057-1067 SEPTEMBER 2021

TB or not TB? Definitive determination of species within the *Mycobacterium tuberculosis* complex in unprocessed sputum from adults with presumed multidrug-resistant tuberculosis

Peter M. Mbelele^{1,2}, Elingarami Sauli², Emmanuel A. Mpolya², Sagal Y. Mohamed³, Kennedy K. Addo⁴, Sayoki G. Mfinanga^{5,6}, Scott K. Heysell³ and Stellah Mpagama^{1,2}

1 Kibong'oto Infectious Diseases Hospital, Kilimanjaro, Tanzania

2 Department of Global Health and Biomedical Sciences, School of Life Sciences and Bioengineering, Nelson Mandela African

Institution of Science and Technology, Arusha, Tanzania

3 Division of Infectious Diseases and International Health, University of Virginia, Charlottesville, VA, USA

4 Department of Bacteriology, Noguchi Memorial Institute for Medical Research, University of Ghana, Accra, Ghana

5 National Institute for Medical Research, Muhimbili Center, Dar es salaam, Tanzania

6 Muhimbili University of Health and Allied Sciences, Dar es salaam, Tanzania

Abstract

OBJECTIVES Differences among Mycobacterium tuberculosis complex (MTC) species may predict drug resistance or treatment success. Thus, we optimised and deployed the genotype MTBC assay (gMTBC) to identify MTC to the species level, and then performed comparative genotypic drugsusceptibility testing to anti-tuberculosis drugs from direct sputum of patients with presumed multidrug-resistant tuberculosis (MDR-TB) by the MTBDRplus/sl reference method. METHODS Patients with positive Xpert[®] MTB/RIF (Xpert) results were consented to provide earlymorning-sputum for testing by the gMTBC and the reference MTBDRplus/sl. Chi-square or Fisher's exact test compared proportions. Modified Poisson regression modelled detection of MTC by gMTBC. RESULTS Among 73 patients, 53 (73%) were male and had a mean age of 43 (95% CI; 40-45) years. In total, 34 (47%), 36 (49%) and 38 (55%) had positive gMTBC, culture and MTBDR respectively. Forty patients (55%) had low quantity MTC by Xpert, including 31 (78%) with a negative culture. gMTBC was more likely to be positive in patients with chest cavity 4.18 (1.31-13.32, P = 0.016), high-quantity MTC by Xpert 3.03 (1.35–6.82, P = 0.007) and sputum smear positivity 1.93 (1.19–3.14, P = 0.008). The accuracy of gMTBC in detecting MTC was 95% (95%) CI; 86-98; $\kappa = 0.89$) compared to MTBDRplus/sl. All M. tuberculosis/canettii identified by gMTB were susceptible to fluoroquinolone and aminoglycosides/capreomycin. CONCLUSIONS The concordance between the gMTBC assay and MTBDRplus/sl in detecting MTC was high but lagged behind the yield of Xpert MTB/RIF. All M. tuberculosis/canettii were susceptible

keywords genotype MTBC, M. tuberculosis complex species, multidrug-resistant, sputum

Sustainable Development Goal: Good health and well-being

to fluoroquinolones, a core drug in MDR-TB treatment regimens.

Introduction

Tuberculosis (TB) is a treatable disease yet is routinely the leading cause of mortality and morbidity from a single infectious pathogen worldwide [1]. TB in human and animals is caused by the species of the *Mycobacterium tuberculosis* complex (MTC), including *M. tuberculosis*, *M. africanum*, *M. bovis*, *M. caprae*, *M. microti*, *M. pinnipedii*, *M. canettii* and *M. mungi* [2]. Previous comparative genomic studies for MTC have described important variations among species in the complex with differences in host association, drug resistance, virulence and epitope diversity [3]. Unfortunately, these studies are predominantly done in countries where TB burden is low but have high capacity for genomic testing and bioinformatics analyses [4]. This may partly limit the understanding of the global distribution of MTC species.

Tanzania is among the top 30 countries with high burden of TB in the world, accounting for 253 cases per 100 000 population incident rate, but a low rifampicin

and/or multidrug-resistant (R/MDR)-TB burden [5]. For instance, in 2018, an estimated 3.4% of new and 18% of previously treated patients for drug-sensitive (DS)-TB had R/MDR-TB globally compared to 1.1% and 3.9% respectively in Tanzania [5]. Of the MTC species, M. tuberculosis is the commonest causative agent for TB in humans globally. The M. africanum is historically predominant in West African countries like Ghana and Senegal [6]. In the past 20 years, one cross-sectional study reported M. bovis as the causative pathogen in 10% of patients with DS-TB in Tanzania [7], but there have been few or no further data, particularly in patients with R/ MDR-TB. Unlike most strains of M. tuberculosis, lack of the pyrazinamidase enzyme in M. bovis leads to natural resistance to pyrazinamide, a key sterilising anti-TB drug for treating patients with DS- and R/MDR-TB [8]. Likewise, a study from West Africa showed that isoniazid resistance due to katG S315T mutation is more frequent in M. tuberculosis than M. africanum [6]. Furthermore, independent of drug resistance, certain species of MTC appear to respond more readily to anti-TB treatment than others [9]. Even with the alarming global trend of drugresistant TB, it is scarcely known whether M. tuberculosis, or M. bovis or M. africanum causes R/MDR-TB. This information is not only important for tracking MTC species circulating in the region, but also in guiding clinical decisions in patients with R/MDR-TB.

The current molecular methods like Xpert® MTB/RIF (Cepheid, USA) and the line probe assays (LPA, Hain LifeScience, Germany) such as the genotype MTBDRplus and MTBDRsl have been endorsed by the World Health Organization (WHO) for detecting MTC and susceptibility to rifampicin, isoniazid, fluoroquinolone and secondline injectable drugs like aminoglycosides/capreomycin [10]. While WHO and a recent clinical trial have recommended use of the MTBDRplus and MTBDRsl in detecting MTC and drug resistance in direct sputum samples collected from patients with R/MDR-TB, they do not detect MTC to the species level [11]. The genotype MTBC assay (gMTBC) is a commercial platform (Hain LifeScience, Germany) that can address this gap, but it is only validated for testing cultured MTC isolates [12,13]. Based on whether mutation(s) are present or absent in 23 rRNA, gyrB and RD1 genes, the gMTBC allows differential identification of the MTC to M. tuberculosis/ M. canettii, M. africanum, M. microti, M. bovis subsp. bovis, M. bovis subsp. caprae and M. bovis BCG only [14]. Similarly, PCR amplification of regions of differences (RDs), spoligotyping or targeted DNA sequencing (i.e. gyrB and hsp65) and whole genome sequencing not only require MTC culture, but also are expensive and largely depend upon extensive bioinformatic expertise

[4,15]. Nevertheless, MTC culture is laborious work and it delays results for up to 8 weeks, yet is prone to contamination in up to 15% of Tanzanian patients [16], and can miss the non-culturable strains [17], and thereby compromise patient care.

While building from the same principle of genotype MTBDRplus and MTBDRsl [11], we optimised and deployed the gMTBC assay to identify MTC to the species level from direct sputum samples of patients with presumed MDR-TB in Tanzania. We also performed comparative susceptibility testing of MTC species identified to rifampicin and isoniazid by using the genotype MTBDRplus, and to fluoroquinolones and second-line injectable drugs by using the genotype MTBDRsl.

Materials and methods

Study settings

This study was conducted from September 2018 through March 2019 at the Kibong'oto Infectious Diseases Hospital (KIDH) in Siha District, Kilimanjaro in northern Tanzania. KIDH is a public hospital with a bed capacity of 320. Each year, the hospital provides services to over 200 and 600 patients with drug-resistant and drug-sensitive TB respectively. It is the National Centre of excellence for clinical management of drug-resistant TB in the country. Patients are usually referred to KIDH if they have R/ MDR-TB with or without comorbidity like diabetes or have confirmed pre/XDR-TB or those living far from R/ MDR-TB ambulatory centres for advanced clinical management. The hospital mycobacteriology laboratory has capacity for testing MTC by smear microscopy for acidfast bacilli (AFB), Lowenstein-Jensen (LJ) solid culture media, Xpert[®] MTB/RIF assay, the line probe assays including the genotype MTBDRplus and MTBDRsl [10], and the gMTBC kits which were added for this study.

Study design and patients

This was a cross-sectional design conducted among patients with presumptive MDR-TB. Presumed MDR-TB was referred to patients who presented with symptoms or signs and risk factors suggestive of R/MDR-TB including prior history of treatment for drug-sensitive TB and HIV with or without rifampicin resistance [18]. The study was approved by the National Ethics Health Review Committee at the National Institute for Medical Research (NIMR) in Tanzania (NIMR/HQ/R.8a/Vol. IX/2662). Permission to conduct the study was granted by KIDH authorities. Patients were included in the study if were aged ≥18 years old but presenting with clinical features

of TB. Prior to any study procedure, they signed a written or oral, witnessed informed consent. A copy of the written consent is available for review by the Editor-in-Chief of this journal on request. Moribund patients and those unable to expectorate quality sputum were excluded. Eligible patients provided early-morningsputum samples for laboratory procedures. Also, socialdemographic and clinical characteristics such as age, gender, duration of sickness at enrolment, HIV status, previous history of TB treatment and working in mining, body weight in kilograms and height in metres were collected from patients. A chest radiograph was also taken to ascertain the presence or absence of lung disease including a cavitary lesion. All patients received anti-TB medications and those living with HIV received anti-retroviral according to the 2018 treatment guidelines in Tanzania.

Sample size estimation and sampling procedure

The minimum sample size for this cross-section study was estimated using 4% prevalence of R/MDR-TB among patients previously treated for DS-TB in Tanzania, with a two-sided type I error of 5%. Assuming that 90% of patients would have agreed to participate in this study, a minimum of 66 patients with presumed R/MDR-TB were required. These patients were conveniently recruited to participate in the study.

Culture, smear microscopy and Xpert® MTB/RIF assay

Each participant provided one early-morning-sputum samples for all laboratory tests. Sputum was decontaminated and cultured in two Lowenstein-Jensen (LJ) slants in according to the Clinical and Laboratory Standard Institute [19] and previous publications [20,21]. Decontaminated sputum was smeared and examined for AFB using light-emitting diode fluorescence microscopy as instructed [22]. All procedures and results interpretations of Xpert[®] MTB/RIF assay were performed as previously described [20,23].

Line probe assays (LPA)

DNA was extracted from decontaminated sputa sediments using the GenoLyse[®] kit in according to manufacturer instructions. The DNA was stored at -20° C until amplification, hybridisation and detection by the genotype MTBDRplus, MTBDRsl and gMTBC kits. Except for gMTBC assay, a multiplex amplification with biotinylated primers was performed on TC 4000 thermal cycler as per prior studies and manufacturer instructions of the genotype MTBDR*plus* and MTBDR*sl* assays. Unlike the previous assay's validation studies and manufacturer instructions for testing clinical isolates [13], the cycling conditions for the genotype MTBC assay were modified to allow testing using direct sputum samples. Briefly, one cycle at 95°C for 15 min, followed by 20 cycles at 95°C for 30 s and 58°C for 2 min in the first stage. In the second stage, this was followed by 30 cycles at 95°C for 25 s, 53°C for 40 s and 70°C for 40 s before a single extension cycle at 70°C for 8min. A M. tuberculosis H37RV DNA reference strain and sterile molecular grade water were run together with DNA extractants, as positive and negative control for MTC respectively. As instructed by manufacturer and a recent reported by Ahmed et al. genotype Mycobacterium CM VER 2.0 kit (Hain LifeScience, Nehren, Germany) was used to test direct sputa samples and cultured isolates that had negative MTC by the MTBDRplus/sl and gMTBC [24]. Amplicons were finally held at 4°C until the DNA stripbased hybridisation and downstream detection steps on twin-incubator, and results were interpreted according to manufacturer instructions and previous publication [25]. Like in Xpert[®] MTB/RIF assay, LPA results are usually available between 24-48 h.

Data management and statistical analysis

A clinical report form was used to collect patient's sociodemographic and clinical data. Bacterial load was categorised as low (sum of low and very low MTC) if MTC was detected at the quantification cycle (Cq) of 23 to >28, and high (sum of high and medium MTC) if MTC was detected at the Cq of < 23 by the Xpert[®] MTB/RIF assay [23]. Chi-square or Fisher's exact test compared proportions. Accordingly, continuous variables were either reported using mean and a 95% confidence interval (CI) or median and their 25th and 75th interguartile range (IQR). Sensitivity, specificity, positive and negative predictive values and accuracy of the genotype MTBC in detecting MTC were calculated as previous [26], using the genotype MTBDRplus/sl, or smear microscopy and LJ culture reference methods. An accuracy was computed as the ratio of the total number of correct results by both test over the total number of patients. The weighted Cohen's kappa (κ) statistic measured the level of agreement between genotype MTBC assays with other tests [27]. Modified Poisson regression was used to model detection of MTC by gMTBC and was adjusted against gender, age, presence of cavitary disease, HIV status, prior history of anti-TB exposure and working in mining with radiological features suggestive of silicosis, smear results for AFB and bacterial load. A 95% confidence interval of test performance was included. A P value

<0.05 was considered significant. All analyses and visualisation were performed in R programming language, version 4.0.2 (http://www.R-project.org).

Results

Socio-demographic and clinical characteristics of patients

Among 73 patients, 73% were male and had an overall mean age of 43 (95% CI: 40-45) years. The sociodemographics and clinical characteristics of patients are presented in Table 1. In total, 39 (53%) had cavity on chest radiograph. Of 39 patients with cavity, 29 (74%) had high MTC quantity compared to 4 (12%) of 34 patients without cavity on chest radiograph (Figure 1a, P < 0.001). Moreover, 27 (37%) were people living with HIV/AIDS (PLWHA) and had a mean CD4 count of 244 (95% CI 184–304) cells/µl. Among PLWHA, 17 (63%) had low MTC quantity compared to 23 (50%) of 46 patients without HIV infection (P = 0.381). Also, 43 (59%) patients had prior history of exposure to first line anti-TB drugs compared to 30 (41%) of anti-TB naive (P = 0.030). However, 20 (67%) of anti-TB naïve patients had low MTC quantity compared to 20 (47%) of those who had prior treatment (P = 0.084), but there was no difference in proportion of culture positivity between patients who had prior treatment and anti-TB naïve [51% (22/43) vs. 47% (14/30), P = 0.887]. Distribution of test positivity/negativity among patients with high and low MTC detected by Xpert MTB/RIF is displayed in Figure 1.

Detection of MTC and non-tuberculous mycobacteria to the species level

Among 73 patients with positive Xpert[®] MTB/RIF assay, 40 (55%) had low MTC quantity. Of these 73 patients, 26 (36%) and 36 (49%) had positive acid-fast bacilli (AFB) on smear microscopy and LJ culture respectively. Also, 38 (52%) and 34 (47%) had MTC detected by genotype MTBDRplus/sl and gMTBC assays, respectively (Figure 2). Overall, 33 (45%) of patient's specimens had positive Xpert® MTB/RIF but negative culture and no MTC detection by genotype MTBDRplus/sl and gMTBC assays. In total, 31 (78%) of 40 patients with low MTC quantity had negative culture, and any of the genotype MTBDRplus/sl and gMTBC compared to 1 (3%) of 33 patients with high MTC quantity (Figure 1d, P < 0.001). In multivariate Poisson regression model, patients with chest cavity, high MTC quantity and smear positive were 4.18 (1.31–13.32, P = 0.016), 3.03 (1.35–6.82, P = 0.007) and 1.93 (1.19–3.14, P = 0.008) times more

Table I Socio-demographic and clinical characteristics (N = 73)

Patient's characteristics	Results		
Male gender, <i>n</i> (%)	53 (73)		
Mean age (95% CI) in years	43 (40-45)		
Age groups in years			
18–30	10 (14)		
31–50	47 (64)		
Above 50	16 (22)		
Median (IQR) days of sickness to enrolment	92 (58-224)		
Cough, n (%)	73 (100)		
Difficulty in breathing, n (%)	53 (73)		
Chest pain, n (%)	58 (79)		
Fever and excessive night sweet, n (%)	54 (74)		
Weight loss of 3 Kg per month, n (%)	58 (79)		
Mean BMI (95% CI) in Kg/m ²	18 (17-20)		
Prior history of worked in mining, n (%)	27 (37)		
HIV positive, n (%)	27 (37)		
Mean (95% CI) CD4 count ($n = 27$)	244 (184-		
	304)		
Prior history of TB treatment, n (%)	43 (59)		
Had cavity on chest radiograph, n (%)	39 (53)		
Low MTC detected by Xpert [®] MTB/RIF, <i>n</i> (%)	40 (55)		

BMI, body mass index; gMTBC, genotype MTBC v1.x; IQR, Interquartile range; MTBDRplus/sl, genotype MTBDRplus or MTBDRsl; MTC, M. tuberculosis complex.

Patients with low MTC included those who had low (5/40) and very low (35/40) MTC at the Xpert® MT/RIF assay cycle quantification (Cq) of 23-28 and >28 respectively. MTC was categorised as high (33/73) if it was detected at the Cq of <23.

likely to have detectable MTC from direct sputa by the gMTBC compared to patients without chest cavity, low MTC and smear negative (Table 2).

Of 36 patients with positive culture, 34 (94%) were identified as MTC and further speciated to M. tuberculosis/canettii by the gMTBC test on both cultured isolates and direct sputa. Neither M. bovis nor M. africanum was detected. The genotype MTBDRplus/sl and gMTBC assays were negative in the remaining 2 (6%) of 36 positive cultures (Figure 2, Figure 1d) and were identified using the genotype Mycobacterium CM VER 2.0 kit as non-tuberculous mycobacteria with one being M. intracellulare and the other one being M. kansasii.

Performance of the genotype MTBC assay in patient's sputum samples

Compared to genotype MTBDRplus/sl, the accuracy of gMTBC in detecting MTC from direct sputa samples was 95% (95% CI; 86–98; Cohen's kappa (κ) = 0.89). Sensitivity, specificity, positive predictive value (PPV), negative predictive (NPV) value and accuracy of the genotype



Figure 1 Distribution of bacterial load measured by Xpert[®] MTB/RIF among patients with chest cavity and HIV and test positivity/ negativity. Bacterial load was quantified as low and high at the Xpert[®] MTB/RIF quantification cycle of 23 to >28, and <23, respectively. Patients were tested using culture and line probe assay (LPA) including the genotype MTBC (gMTBC) and genotype MTBDRplus/sl (MTBDR).

MTBC assay in detecting MTC from direct patient's sputa were summarised in Table 3. There was a strong concordance between the gMTBC assay and both LJ culture and MTBDR*plus/sl* ($\kappa \ge 89\%$), but was moderate when compared to smear microscopy ($\kappa = 61\%$, Table 3). Except for 2 patients with non-tuberculous mycobacteria, the gMTBC detected MTC in 33 patients with positive culture at an accuracy of 95% (95% CI; 86–98) compared to culture (Figure 2 and Table 3). Moreover, gMTBC detected 7 patients more than smear microscopy (Figure 2).

Susceptibility profile of MTC to anti-TB drugs

Among 38 patients with MTC, 28 (74%) had rifampicin resistance detected by both Xpert[®] MTB/RIF and genotype MTBDRplus assays. Eight (29%) of 28 patients with rifampicin-resistant MTC strains had isoniazid resistance, suggestive of MDR-TB. One patient had isoniazid resistance in MTC but rifampicin was susceptible. All 38 patients were susceptible to fluoroquinolones and injectable aminoglycosides/peptides as shown by genotype MTBDRsl.

Discussion

This study showed strong concordance between the genotype MTBC and MTBDRplus/sl as well as the LJ culture, in detecting MTC from sputum samples of patients with presumed R/MDR-TB. The assay achieved a $\geq 97\%$ sensitivity compared to LJ culture and genotype MTBDRplus/ sl which is higher than 93% sensitivity reported by Somoskovi *et al.*, from patients with culture and smearpositive sputa [13]. Molecular methods including the gMTBC and sequencing technologies have been optimised for testing high burden samples such as cultured isolates



Figure 2 Detection of *M. tuberculosis* complex (MTC) from patients presumed for multidrug-resistant (MDR)-TB. MTC are the causative agents for tuberculosis (TB), and were identified using the line probe assays including the genotype MTBC (gMTBC), and any of the genotype MTBDRplus or MTBDRsl assays (MTBDR), and were further speciated to *M. tuberculosis/canettii* by gMTBC. Non-tuberculous mycobacteria (NTM) were identified to *M. intracellulare and M. kansansii* by using the genotype Mycobacterium CM VER 2.0 kit.

[13,25]. Despite that the genotype MTBDRplus and MTBDRsl assays have recently been tested on unprocessed sputa [11], this is to our knowledge the first report to compare the performance with the gMTBC.

Patients with cavitary disease had high MTC compared to those without chest cavity. Presence of cavity can also result in an local environment of high bacterial replication within the cavity which can later distribute into the airway, consistent with a higher sputum bacterial burden [28,29]. Patients with chest cavity, HIV negative, smear positive, and high bacterial load measured by Xpert® MTB/RIF harboured detectable MTC from direct sputa by the gMTBC compared to patients without chest cavity, low MTC and smear negative. Sensitivity improves in people with conditions of known higher bacterial burden and expectedly decreases in paucibacillary conditions particularly in patients with atypical features of TB on chest radiograph [30,31]. Moreover, the high yield of gMTBC compared to smear microscopy supports its use during screening of participants in clinical trials [11].

Surprisingly in this study, approximately 45% of patients with a positive Xpert[®]MTB/RIF result for MTC had a negative LPA and LJ culture, as compared to only 20% from other studies [32]. Under expected circumstances, culture should be more sensitive than Xpert[®] MTB/RIF assay and LPA [33]. A similar discordance was

recently reported in a clinical trial in which 34% of paucibacillary samples re-tested negative for active TB by Xpert[®] MTB/RIF assay [34,35]. Similar to our findings, a paucibacillary state was among the potential sources of discrepancy [34,35]. Despite the similar culture positivity rate between patients with prior history of treatment with first-line anti-TB drugs and those who were anti-TB drug naïve, approximately half of those with prior treatment had low detectable MTC by Xpert® MTB/RIF assay, which may have represented a false-positive signal [31]. Moreover, in patients with prior history of TB treatment, detecting low level of MTC DNA by Xpert® MTB/RIF or the more sensitive Xpert® MTB/RIF Ultra may represent a subpopulation without active infection but rather chronic lung disease, and residual cavitary changes and/ or bronchiectasis [36,37]. In such patients, the differential diagnosis should also include non-tuberculous mycobacteria and fungal pathogens such as Aspergillus species that are known to colonise those with residual lung disease following otherwise successful TB treatment [37,38]. We and others have reported on the importance of accurately detecting non-tuberculous mycobacteria such as M. intracellulare and M. kansasii among people living with HIV/ AIDS in Tanzania for whom multiple rounds of TB treatment were performed before alternative diagnoses and post-TB chronic lung disease were considered [39,40]. To enhance more accurate TB diagnosis, our findings

	Univariate		Multivariate		
Variable	Crude incident rate ratio (95% CI)	P value	Adjusted incident rate ratio (95% CI)	P value	
Gender					
Female	Ref:		Ref:		
Male	0.79 (0.38-1.61)	0.518	1.10 (0.66–1.82)	0.725	
Age group in years	3				
18-30	Ref:		Ref:		
31-50	0.61 (0.35-1.08)	0.093	1.58 (1.11-2.27)	0.012	
Above 50	0.96 (0.49–1.88)	0.899	0.49 (0.31–0.78)	0.003	
HIV status					
Negative	Ref:		Ref:		
Positive	0.93 (0.46-1.88)	0.838	0.53 (0.36-0.79)	0.002	
Cigarette smoking					
Non-smokers	Ref:		Ref:		
Smokers	1.04 (0.53-2.03)	0.918	0.81 (0.55-1.19)	0.283	
Chest cavity	X Z				
Absent	Ref:		Ref:		
Present	9.52 (2.91-31.13)	< 0.001	4.18 (1.31–13.32)	0.016	
MTC quantity					
Low	Ref:		Ref:		
High	7.03 (2.72–18.16)	< 0.001	3.03 (1.35-6.82)	0.007	
Chest silicosis					
Absent	Ref:		Ref:		
Present	1.05 (0.53-2.11)	0.88	0.96 (0.63-1.47)	0.863	
Prior TB treatment	t				
New	Ref:		Ref:		
Retreatment	1.13 (0.56-2.25)	0.735	0.84 (0.61-1.17)	0.231	
Smear microscopy					
Negative	Ref:		Ref:		
Positive	3.78 (1.84–7.75)	< 0.001	1.93 (1.19-3.14)	0.008	

Table 2 Modified Poisson regression modelling of MTC detection by gMTBC (N = 73)

gMTBC, genotyope MTBC ver 1.x; HIV, human immunodeficiency virus; MTC, M. tuberculosis complex measured by Xpert[®] MTB/ RIF; TB, tuberculosis.

warrant consideration of modifying the current Xpert[®] MTB/RIF testing algorithm to include an assessment of active/replicating bacillary load with tools such as the Molecular Bacterial Load assay, which detects and quantifies the rRNA of viable M. tuberculosis within 24 h, as well as other conventional techniques such as culturing the M. tuberculosis in liquid and solid media for those with prior TB treatment, and expanding access to nonsputum assays such as urine TB-LAM for people living with HIV/AIDS, and suspicion of concomitant extrapulmonary TB [41-43]. In doing so, patients with active TB can be treated with anti-TB drugs, while those with nontuberculous mycobacteria or other manifestations of post-TB lung disease can be effectively triaged for chronic lung disease pharmacotherapy and pulmonary rehabilitation [44]. These approaches support recent observations such as those by Costantini et al (2020) who reported a favourable outcome after withholding anti-TB medication

in a patient with positive Xpert[®] MTB/RIF but negative culture results [45].

The M. tuberculosis/canettii was the only member of MTC identified as the pathogen among patients with presumptive MDR-TB. Because of the close genetic relatedness, the gMTBC assay cannot fully differentiate M. tuberculosis from M. canettii, even with the recent 894 diverse genomes of M. canettii and major phylogenetic groups of MTC that were typed to guide the new revised interpretation of the assay's result [46]. The M. tuberculosis/cannetii predominance in this study is similar to a prevalence survey in sub-Sahara Africa, that at least 95% of TB was due to M. tuberculosis and 2% due to M. africanum [25]. Our findings confirm prior epidemiological studies that had suggested the geographic distribution of M. africanum shows it to be an important cause of TB in West Africa, but that it is rarely detected in East African countries, including Tanzania [47]. In

Reference methods	Genotype	Genotype MTBC assay, n % (95% confidence interval)								
	Positive	Negative	Sensitivity	Specificity	PPV	NPV	Accuracy	Cohen's kappa		
Smear Microscopy										
Positive	23	3	68 (49-83)	92 (78-98)	88 (69–97)	77 (62-87)	81 (70-89)	61 (43–79)		
Negative	11	36								
LJ culture										
Positive	33	3	97 (83-100)	92 (78–98)	92 (66–99)	97 (84-100)	95 (86-98)	89 (79–99)		
Negative	1	36								
Genotype MTBDR	olus/sl									
Positive	34	4	100 (87-100)	90 (75–97)	89 (74–97)	100 (88-100)	95 (86-98)	89 (79–99)		
Negative	0	35								

Table 3 Performance of the genotype MTBC assay in detecting M. tuberculosis from direct sputum samples

LJ, Lowenstein-Jensen solid media; MTBDRplus/sl, genotype MTBDRplus or MTBDRsl assays for detecting *M. tuberculosis* complex and drug resistance; NPV, negative predictive value; PPV, positive predictive value.

these settings, this species dominance underscores the treatment, prevention and control strategies for *M. tuber-culosis*.

Isoniazid resistance in patients with rifampicin resistant was uncommon, accounting for about 29% compared to over 64% in South Africa [48] and 78% in China [49]. In the present study, isoniazid monoresistance was categorised from direct sputum samples which was negative in a large number of patients [30], as compared to MTC isolates in two studies [48,49]. However, this is contrary to a global genomic analysis of over 5000 M. tuberculosis strains by Manson et al., who documented that mutations conferring resistance to isoniazid evolve first before rifampicin resistance across all lineages, geographic regions and time periods [50], and therefore our findings likely represent the incomplete coverage of isoniazid resistance targets on the genotype MTBDRplus assay [51]. We were not able to perform conventional culture-based isoniazid susceptibility testing to interrogate this observation. Favourably in the present study, all M. tuberculosis/canettii identified were susceptible to fluoroquinolone and aminoglycosides/capreomycin, similar to findings by Mpagama et al., among patients with MDR-TB in Tanzania [52] and in studies from Kenya [53]. Our findings continue to support the empirical use of fluoroquinolones in the treatment of R/MDR-TB in Tanzania while awaiting treatment modification from culture-based susceptibility testing.

This study has limitations. The relatively small sample size does not allow for generalisability beyond the referral patterns to the hospital of study, but do provide a confirmatory snapshot of circulating MTC species in patients with MDR-TB in Tanzania to inform further the current approach to optimal diagnosis and clinical management. However, the main comparison was made in the performance characteristics in gMTBC from unprocessed sputa that did not require a diversity of samples. Furthermore, we included patients with previous history of treatment or those with rifampicin-resistant TB rather than anti-TB naïve, which could have contributed to discrepancies among tests as previously outlined [54]. Nonetheless, these patients represent an important population at risk for acquiring drug resistance or other circulating strains given their exposure to the healthcare setting and as evidenced by our findings of nontuberculous mycobacteria.

In conclusion, the concordance between the gMTBC and MTBDRplus/sl in detecting MTC from direct sputum was considerably high but lagged behind the yield of Xpert MTB/RIF. High bacterial loads such as those found in patients with chest cavity and sputum smear positive were independent predictors of MTC detection. Modification by integrating biomarkers for active TB, such as the TB molecular bacterial load assay, urine LAM and selective use of mycobacterial culture for nontuberculous mycobacteria into Xpert® MTB/RIF testing algorithm for those with prior TB treatment can be pivotal not only in addressing potentially false-positive Xpert[®] MTB/RIF results, but also in differentiating active TB from the sequelae of post-TB lung disease. Importantly, all M. tuberculosis/canettii detected were susceptible to fluoroquinolones, a core drug in MDR-TB treatment regimens.

Acknowledgements

The authors acknowledge to Dr Bibie Saidi for administering informed consent and questionnaire; Ms Batuli Mono, Ms Taji Mnzava, Ms Rovina Kiwia, and Mr Joseph Kachalla for administering informed consent and

sputum sample collection from participants; and Ms Sarapia Mallya and Mr Hussein Mfugale for handling sputum samples and their assistance during performing mycobacteriology procedures. Authors also acknowledge Dr. Katharina Krepple for proof reading the article. In addition, authors thank the KIDH administrative authorities for granting permission to conduct this study.

Funding information

The authors acknowledge funding support from the DEL-TAS Africa Initiative (Afrique One-ASPIRE/DEL-15-008). Afrique One-ASPIRE is funded by a consortium of donors including the African Academy of Sciences, Alliance for Accelerating Excellence in Science in Africa, the New Partnership for Africa's Development Planning and Coordinating Agency, the Wellcome Trust (107753/A/15/ Z), and the UK Government. Likewise, PMM was supported by the EDCTP2 programme supported by the European Union (grant number TMA2016SF - SF1463 — REMODEL-TZ, a senior fellowship for SGM). Moreover, R01 AI137080 supported to SKH. However, these funding bodies had no role in the idea conception and design, collection, analysis, and interpretation of data and in writing of manuscript.

Data availability statement

The data used in this manuscript are available from the corresponding author upon request and permission from Kibong'oto Infectious Diseases Hospital (KIDH).

References

- 1. World Health Organization. *Global tuberculosis report* 2018. Switzerland: Geneva, Switzerland, 2018.
- Riojas MA, Mcgough KJ, Rider-riojas CJ, Rastogi N, Hazbón MH. Phylogenomic analysis of the species of the Mycobacterium tuberculosis complex demonstrates that Mycobacterium africanum, Mycobacterium bovis, Mycobacterium caprae, Mycobacterium microti and Mycobacterium pinnipedii are later heterotypic synonyms of My. Int J Syst Evol Mycrobiology 2018: 68: 324–332.
- 3. Jia X, Yang L, Dong M *et al*. The bioinformatics analysis of comparative genomics of *Mycobacterium tuberculosis* complex (MTBC) provides insight into dissimilarities between intraspecific groups differing in host association, virulence, and epitope diversity. *Front Cell Infect Microbiol* 2017: 7: 1–14.
- 4. Colman RE, Mace A, Seifert M *et al*. Whole-genome and targeted sequencing of drug-resistant *Mycobacterium tuber-culosis* on the iSeq100 and MiSeq: a performance, ease-of-

use, and cost evaluation. *PLoS Medicine* 2019: 16: e1002794.

- World Health Organization. *Global tuberculosis report* 2019. World Health Organization: Geneva, Switzerland, 2019.
- Otchere ID, Asante-Poku A, Osei-Wusu S *et al*. Detection and characterization of drug-resistant conferring genes in *Mycobacterium tuberculosis* complex strains: a prospective study in two distant regions of Ghana. *Tuberculosis* 2016: 99: 147–154.
- Mfinanga SG, Morkve O, Kazwala RR et al. Mycobacterial adenitis: role of Mycobacterium bovis, non-tuberculous mycobacteria, HIV infection, and risk factors in Arusha, Tanzania. East Afr Med J 2004: 81: 171–178.
- Valle MB, Torres-gonzález P, Cervera-hernández ME. Trends of *Mycobacterium bovis* isolation and first-line antituberculosis drug susceptibility profile : a fifteen-year laboratory-based surveillance. 2015, 1–11.
- Scott C, Cavanaugh JS, Silk BJ et al. Comparison of sputum-culture conversion for Mycobacterium bovis and M. tuberculosis. Emerg Infect Dis 2017: 23: 456–462.
- Mbelele PM, Mohamed SIY, Sauli E *et al*. Meta-narrative review of molecular methods for diagnosis and monitoring of multidrug-resistant tuberculosis treatment in adults. *Int J Mycobacteriology* 2018: 7: 299–309.
- Alipanah N, Shete PB, Nguyen H *et al.* Feasibility of direct sputum molecular testing for drug resistance as part of tuberculosis clinical trials eligibility screening. *Diagnostics* 2019: 9: 1–5.
- Safianowska A, Walkiewicz R, Nejman-Gryz P, Chazan R, Grubek-Jaworska H. Diagnostic utility of the molecular assay GenoType MTBC (HAIN Lifescience, Germany) for identification of tuberculous mycobacteria. *Pneumonol Alergol Pol* 2009: 77: 517–520.
- Somoskovi A, Dormandy J, Rivenburg J, Pedrosa M, McBride M, Salfinger M. Direct comparison of the Geno-Type MTBC and genomic deletion assays in terms of ability to distinguish between members of the *Mycobacterium tuberculosis* complex in clinical isolates and in clinical specimens. J Clin Microbiol 2008: 46: 1854–1857.
- 14. LifeSciences H. Geno type MTBC Ver 1.X: instructions for use, IFU-301-09. GmbH: Nehren, Germany, 2012.
- 15. Shea J, Halse TA, Lapierre P *et al.* Comprehensive wholegenome sequencing and reporting of drug resistance profiles on clinical cases of *Mycobacterium tuberculosis* in New York State. *J Clin Microbiol* 2017: 55: 1871–1882.
- Hoza AS, Mfinanga SGM, König B. Anti-TB drug resistance in Tanga, Tanzania: a cross sectional facility-base prevalence among pulmonary TB patients. *Asian Pac J Trop Med* 2015: 8: 907–913.
- Barac A, Karimzadeh H, Mahya E, Mohammad P, Rahimi T. Laboratory cross - contamination of *Mycobacterium tuberculosis*: a systematic review and meta - analysis. *Lung* 2019; 197: 651–661.
- 18. World Health Organization. Companion handbook to the WHO guidelines for the programmatic management of

drug-resistant tuberculosis. World Health Organization: Geneva, Switzerland, 2014.

- Clinical and Laboratory Standards Institute (CLSI). Laboratory detection and identification of mycobacteria; approved standard. CLSI document M48-A. Clinical and Laboratory Standards Institute: Wayne, PA, 2008.
- Mbelele PM, Aboud S, Mpagama SG, Matee MI. Improved performance of Xpert MTB/RIF assay on sputum sediment samples obtained from presumptive pulmonary tuberculosis cases at Kibong'oto infectious diseases hospital in Tanzania. BMC Infect Dis 2017: 17: 1–7.
- Tripathi K, Tripathi PC, Nema S, Shrivastava AK, Dwiwedi K. Modified petroff 's method : an excellent simplified decontamination technique in comparison with petroff 's method. *Int J Recent Trends Sci Technol* 2014: 10: 461–464.
- Lumb R, Van DA, Bastian I, Fitz-Gerald M. Laboratory diagnosis of tuberculosis by sputum microscopy. 2013: (Available from: http://www.stoptb.org/wg/gli/assets/docume nts/tb%20microscopy%20handbook_final.pdf.) [19 Jun 2021]
- Ssengooba W, Respeito D, Mambuque E et al. Do Xpert MTB/RIF cycle threshold values provide information about patient delays for tuberculosis diagnosis ? PLoS One 2016: 302: 1–10.
- Ahmed ST. USING GenoType Mycobacterium CM/assay for rapid detection and differentiation of non-tuberculous from sputum specimen. *Biochem Cell Arch* 2020: 20: 5973– 5977.
- Addo KK, Addo SO, Mensah GI, Mosi L, Bonsu FA. Genotyping and drug susceptibility testing of mycobacterial isolates from population-based tuberculosis prevalence survey in Ghana. BMC Infect Dis 2017: 17: 1–7.
- 26. Mishra H, Reeve BWP, Palmer Z et al. Xpert MTB/RIF Ultra and Xpert MTB/RIF for diagnosis of tuberculosis in an HIV-endemic setting with a high burden of previous tuberculosis: a two-cohort diagnostic accuracy study. Lancet Respir Med 2020: 8: 368–382.
- 27. McHugh ML. Lessons in biostatistics interrater reliability: the kappa statistic. *Biochem Medica* 2012: 22: 276–282. (Available from: https://hrcak.srce.hr/89395.)
- Ong CWM, Elkington PT, Friedland JS. Tuberculosis, pulmonary cavitation, and matrix metalloproteinases. Am J Respir Crit Care Med 2014: 190: 9–18.
- Murthy SE, Chatterjee F, Crook A *et al.* Pretreatment chest x-ray severity and its relation to bacterial burden in smear positive pulmonary tuberculosis. *BMC Med* 2018: 16: 1– 11.
- Theron G, Venter R, Smith L *et al.* False-positive Xpert MTB/RIF results in retested patients with previous tuberculosis: frequency, profile, and prospective clinical outcomes. J Clin Microbiol 2018: 56: e1696–e1717.
- Theron G, Venter R, Calligaro G *et al.* Xpert MTB/RIF results in patients with previous tuberculosis: can we distinguish true from false positive results? *Clin Infect Dis* 2016: 62: 995–1001.

- 32. Bhardwaj A, Khan S, Kumar A, George L, Mehta A, Radhakrishnan K. Assessing the utility of genexpert MTB/Rif assay in a tertiary care centre in southern India with established microscopy and liquid culture facilities. J Assoc Physicians India 2019: 67: 31–34.
- Rasheed W, Rao NA, Adel H, Baig MS, Adil SO. Diagnostic accuracy of Xpert MTB/RIF in sputum smear-negative pulmonary tuberculosis. *Cureus* 2019: 11: e5391.
- 34. Ngabonziza JCS, Decroo T, Migambi P et al. Prevalence and drivers of false-positive rifampicin-resistant Xpert MTB/ RIF results: a prospective observational study in Rwanda. *Lancet Microbe* 2020: 1: e74–e83.
- Variava E, Moloantoa T, Martinson NA. Discrepancies in Xpert tuberculosis testing. *Lancet Microbe* 2020: 1: e47– e48.
- Allwood BW, Stolbrink M, Baines N et al. Persistent chronic respiratory symptoms despite TB cure is poorly correlated with lung function. Int J Tuberc lung Dis 2021: 25: 262– 270.
- Hsu D, Irfan M, Jabeen K *et al.* Post tuberculosis treatment infectious complications. *Int J Infect Dis* 2020: 92: S41– S45.
- Wanner A, Edwards M, Harries AD *et al.* International research and guidelines on post-tuberculosis chronic lung disorders: a systematic scoping review. *BMJ Glob Heal* 2018: 3: 1–8.
- 39. Mpagama SG, Heysell SK, Ndusilo ND *et al.* Diagnosis and interim treatment outcomes from the first cohort of multidrug-resistant tuberculosis patients in Tanzania. *PLoS One* 2013: 8: e62034.
- Mnyambwa NP, Kim D-J, Ngadaya E et al. Genome sequence of Mycobacterium yongonense RT 955–2015 isolate from a patient misdiagnosed with multidrug-resistant tuberculosis: first clinical detection in Tanzania. Int J Infect Dis 2018: 71: 82–88.
- Honeyborne I, McHugh TD, Phillips PPJ *et al.* Molecular bacterial load assay, a culture-free biomarker for rapid and accurate quantification of sputum *Mycobacterium tuberculosis* bacillary load during treatment. *J Clin Microbiol* 2011: 49: 3905–3911.
- Sabiiti W, Azam K, Farmer ECW *et al.* Tuberculosis bacillary load, an early marker of disease severity: the utility of tuberculosis molecular bacterial load assay. *Thorax* 2020: 1–3.
- 43. Mbelele PM, Mpolya EA, Sauli E *et al.* Mycobactericidal effect of different regimens measured by molecular bacterial load assay among people treated for multidrug-resistant tuberculosis in Tanzania. *J Clin Microbiol* 2021: **59**: e02927–e3020.
- 44. Singh SK, Naaraayan A, Acharya P, Menon B, Bansal V, Jesmajian S. Pulmonary rehabilitation in patients with chronic lung impairment from pulmonary tuberculosis. *Cur*eus 2018: 10: e3664.
- Costantini L, Marando M, Gianella P. Long-term GeneXpert positivity after treatment for pulmonary tuberculosis. *Eur J Case Reports Intern Med* 2020: 7: 1737.

- Loiseau C, Brites D, Moser I *et al*. Revised interpretation of the Hain lifescience genotype MTBC to differentiate mycobacterium canettii and members of the mycobacterium tuberculosis complex. *Antimicrob Agents Chemother* 2019: 63: e159–e219.
- 47. Zumla A, Darko I, Ivy G *et al.* Learning from epidemiological, clinical, and immunological studies on *Mycobacterium africanum* for improving current understanding of host pathogen interactions, and for the development and evaluation of diagnostics, host-directed therapies, and vac. *Int J Infect Dis* 2017: 56: 126–129.
- 48. Pitso L, Potgieter S, Van der Spoel van Dijk A. Prevalence of isoniazid resistance-conferring mutations associated with multidrug-resistant tuberculosis in free state province, South Africa. S Afr Med J 2019: 109: 659–664.
- 49. Liu Z, Zhang M, Wang J *et al*. Longitudinal analysis of prevalence and risk factors of rifampicin-resistant tuberculosis in Zhejiang, China. *Biomed Res Int* 2020: 2020: 1–8.
- 50. Manson AL, Cohen KA, Abeel T et al. Genomic analysis of globally diverse Mycobacterium tuberculosis strains provides

insights into emergence and spread of multidrug resistance. *Nat Genet* 2017: **49**: 395–402.

- 51. Jo KW, Yeo Y, Sung H, Kim MN, Shim TS. Analysis of discrepant results between the Genotype® MTBDRplus assay and an antimicrobial drug susceptibility test for isoniazidresistant tuberculosis. *Respir Med* 2017: 122: 12–17.
- 52. Mpagama SG, Houpt ER, Stroup S et al. Application of quantitative second-line drug susceptibility testing at a multidrug-resistant tuberculosis hospital in Tanzania. BMC Infect Dis 2013: 13: 432.
- 53. Ogari CO, Nyamache AK, Nonoh J, Amukoye E. Prevalence and detection of drug resistant mutations in *Mycobacterium tuberculosis* among drug naïve patients in Nairobi, Kenya. *BMC Infect Dis* 2019: **19**: 279.
- Jaworski JP, Pluta A, Rola-Łuszczak M et al. Interlaboratory comparison of six real-time PCR assays for detection of bovine leukemia virus proviral DNA. J Clin Microbiol 2018: 56: 1–11.

Corresponding Author Peter M. Mbelele, Kibong'oto Infectious Diseases Hospital (KIDH), P.O BOX 12, Siha, Kilimanjaro, Tanzania. E-mail: mbelelepeter@yahoo.com